


Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
Date: \_\_\_\_\_

## circaid juxtafit essentials arm ready-to-wear size chart



circumferences	1	2	3	4	5	6
eG	28 – 32	30 – 34	32 – 37	35 – 40	38 – 44	42 – 48
eF	26 – 30	28 – 32	30 – 35	33 – 38	36 – 41	39 – 45
eE	25 – 28	27 – 30	29 – 32	31 – 34	33 – 36	36 – 40
eD	23 – 26	25 – 28	27 – 30	29 – 32	31 – 34	33 – 37
eC	15 – 17	17 – 19	19 – 21	21 – 23	23 – 25	25 – 28

actual product length	x-short	short	long
lC – G	38	43	48
arm length	x-short	short	long
lC – E	< 22	22 – 25	> 25

Measurements in cm